






STANDARDS | TOOLS
FOR THERMOFORMED UNIDOSE




w w w . s i c o f o r . c o m




GAT OT 4001




GAT OT 4002




GAT OT 4003




GAT OT 4004




GAT OT 4005




GAT OT 4006




GAT OT 4007




GAT OT 4008




GAT OT 4009




GAT OT 4010



GAT OT 4011



SOG OT 8001




PIL OT 4001


RECAPITULATORY TABLE

STANDARDS


REFERENCE	TOOLS	TYPE OF OPENING	DOSING VOLUME	NBER OF LAYING	CUTTING DIMENSION
GAT OT 4001	Jar	Peelable	10 ml	4	72 x 43
GAT OT 4002	Jar	Scored	3 ml	4	72 x 43
GAT OT 4003	Jar	Scored	3 ml	4	72 x 43
GAT OT 4004	Jar	Scored	3 ml	4	72 x 43
GAT OT 4005	Bottle	Scored	3 ml	4	72 x 43
GAT OT 4006	Jar	Scored	3 ml	4	72 x 43
GAT OT 4007	Tube	Peelable	7 ml	4	72 x 43
GAT OT 4008	Jar	Scored	3 ml	4	72 x 43
GAT OT 4009	Jar	Scored	3 ml	4	72 x 43
GAT OT 4010	Tappet	Scored	6 ml	4	72 x 43
GAT OT 4011	Drop	Scored	3 ml	4	72 x 43
SOG OT 8001	Drop	Scored	1,5 ml	8	38 x 40
PIL OT 4001	Dose	Scored	7,5 ml	4	80 x 43




SIC OT 4001




SHI OT 4001



ARK OT 3001




PAT OT 1201




SIC OT 8001




LDB OT 8001




SHI OT 6001




IMA OT 6001




GAT OT 6002




GAT OT 6003



GAT OT 6004



MUG OT 6001



MUG OT 6002

RECAPITULATORY TABLE | STANDARDS

REFERENCE	TOOLS	TYPE OF OPENING	DOSING VOLUME	NBER OF LAYING	CUTTING DIMENSION
SIC OT 4001	Jar	Scored	10 ml	4	72 x 43
SHI OT 4001	Jar	Scored	5 ml	4	80 x 43
ARK OT 3001	Tube	Scored	3 et 5 ml	3 trios	95 x 71
PAT OT 1201	Bubble	Scored	0,4 ml	12	35 x 25
SIC OT 8001	Bubble	Peelable	5 ml	8	38 x 40
LDB OT 8001	Bubble	WITHOUT	WITHOUT	8	WITHOUT
SHI OT 6001	Jar	Scored	5 ml	6	59 x 55
IMA OT 6001	Blister pack	Scored	OBJET	6	88 x 32
GAT OT 6002	Tube	Scored	3 ml	6	WITHOUT
GAT OT 6003	Tube	Scored	5 ml	6	WITHOUT
GAT OT 6004	Tube	Scored	3 ml	6	WITHOUT
MUG OT 6001	Jar	Peelable	3 ml	6 (2x3)	100 x 70
MUG OT 6002	Jar	Peelable	5 ml	6 (2x3)	100 x 70